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a Novartis company

 Splendor
Ciba Vision / Alcon Vision Care distributor

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2. U& A Contact Lenses na Make-Up with ECPs in France & UK. N=200. June-July 2011.

Neke činjenice o make-up-u, proizvodima za skidanje šminke i dr. kozmetičkim proizvodima i kontaktnim lećama

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*Tsukiyama T et al. Influence of eye cosmetics and cleansing products on contact lenses. J Japan Contact Lens Soc 2010; 52:101-7.

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
Figure 1: Contamination of upper lid margin by eye make-up applied to the lower lid margin. (Image courtesy of Alison Ng – Cardiff University)

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
Figure 2: Cosmetic 'glitter' contaminants suspended in pre-corneal tear film (Image courtesy of Alison Ng – Cardiff University)

Često se make-up proizvodi označavaju kao vodootporni, te takva kozmetika sadrži ulja, vosak ili silikonske agense, kako bi imala bolju otpornost na suze i znoj.

Za sada nema objavljenih radova koji ukazuju na efekte kontaminacije make-upom suza, i utjecajem na normalnu funkciju i fiziologiju suznog filma.

Ali često se mogu vidjeti ovakve slike na biomikroskopu – sitne pigmentirane čestice koje plutaju u suznom filmu.

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
Figure 3: Mascara laden rheum (Image courtesy of Alison Ng – Cardiff University)

Zamijećena su relativno velike količine Migracija** kozmetičkih proizvoda za oči već 30 minuta nakon aplikacije make-up proizvoda, a posebno kod onih koji nanose kozmetiku blizu ili na liniju trepavica.

- Studije impliciraju potencijalne promjene u fiziologiji oka, koje utječu na suzni film, te postoje indicije o povezanosti upotrebe Make-up-a sa MGD-om (disfunkcija meibomove žlezde)*
- Utjecaj na nošenje kontaktnih leća nije dovoljno ispitan

* Arnold Silk AC. Eye make-up and the soft contact lens wearer. Ophthal Optician 1982; 2:101-6.

** Goto T et al. Cosmetic product migration onto the ocular surface: Exacerbation of migration after eye drop instillation. Cornea 2010; 29:140-3.

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



Chart 1 For surface treatment, AIR OPTIX® and Night & Day® from CIBA VISION use plasma treatment; ACUVUE OASYS® uses HydroClear® Plus; PureVision® uses Permalon® Surface process.

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| | Change from baseline | | | | |
|-------------------------|----------------------|----------------|---------------------|---------|---------------|
| | Diameter | Sagittal depth | Base curve | Power | Image quality |
| Makeup remover (liquid) | ≤ 0.50 mm (flatter) | ≤ 0.15 mm | ≤ 0.77 mm (flatter) | < 0.25D | Strong impact |
| Makeup remover (wipes) | ≤ 0.10 mm | ≤ 0.07 mm | ≤ 0.13 mm | < 0.25D | Little impact |
| Hand cream | ≤ 0.15 mm | ≤ 0.08 mm | ≤ 0.15 mm | < 0.25D | Little impact |
| Mascara (waterproof) | ≤ 0.10 mm | ≤ 0.07 mm | ≤ 0.12 mm | < 0.25D | Strong impact |
| Mascara (water soluble) | ≤ 0.14 mm | ≤ 0.07 mm | ≤ 0.14 mm | < 0.25D | Strong impact |
| ISO tolerances | ±0.20 mm | ±0.20 mm | ±0.05 mm | ±0.25 D | n/a |

Overall, balafilcon A, enfilcon A and galyfilcon A lenses demonstrated a stronger change in parameters, while plasma surface lotrafilcon A and lotrafilcon B seemed to be impacted the least by the cosmetics.

•The impact of cosmetics on the physical dimension and optical performance of silicone hydrogel contact lenses; Doerte Luensmann, Centre for Contact Lenses Research; University of Waterloo, Canada

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